

**GL4U: Introduction 2024 Lecture 4 of 4** 

Amanda M. Saravia-Butler, Ph.D.

NASA GeneLab Science Lead

Contractor: KBR

Biological & Physical Sciences



## Sequencing Through the Ages




1977

1990



Now






Assembly required

**Assembly** required

Not much assembly required




## Pyrosequencing aka Sequencing by Synthesis



## Illumina Sequencing Workflow


- Illumina is a very successful biotech company specializing in next generation technology that uses the pyrosequencing method
- ~90% of all sequencing worldwide is performed on an Illumina instrument (including GeneLab)
- The Illumina sequencing workflow has the following 3 steps:
  - **➤** Library Construction
  - > Cluster Formation
  - Sequencing




## **DNAseq: Library Preparation**


Step 1: Create DNA fragments (with a means to attach adapters) from the extracted sample DNA

**Tagmentation: Illumina DNA Prep (formerly Nextera Flex)** 






### **Target PCR: Illumina Amplicon**




## **DNAseq: Library Preparation**

### **Step 2: Attach adapters**

**Tagmentation: Illumina DNA Prep** 







## **DNAseq: Library Preparation**

### **Step 3: Amplify libraries**


**Tagmentation: Illumina DNA Prep** 



### Fragmentation/A-tailing: IDT




#### **Target PCR: Illumina Amplicon**




### DNA Library QC: Electrophoresis and Electropherogram

### **Electrophoresis**




### Electropherogram





## **DNAseq: Library QC**

Libraries are evaluated using a bioanalyzer or a tape station to create an electropherogram to assess quality

#### **Good** Library - 1








#### When assessing library quality look for the following:


- NO adapter dimers! Why?
- Library size is consistent with the number of desired sequencing cycles
  - ➤ If you're sequencing at PE 250, what is a good library size? (hint: library size = insert length + adapter length)


Assuming adapters are ~65bp each (130bp) ~600bp that would give an insert length of ~470, allowing a ~30bp overlap between R1 and R2

## **DNAseq: Library QC**

- Assess library quality (bioanalyzer, TapeStation)
  - Adapter dimers
  - > Fragment size
- Determine library quantity (Qubit, qPCR)









% AD	% PF	% AD Reads	
Control	69.54	0.24	
10%	10.87	84.25	
5%	21.39	60.44	
1%	51.88	6.46	



### **rRNA** Contamination

### mRNA makes up only ~2-5% of a total RNA sample





## **RNAseq: Library Preparation**

Figure 1 Ribo-Zero Depleting and Fragmenting RNA



Figure 2 Synthesizing First Strand cDNA




Figure 3 Synthesizing Second Strand cDNA







Figure 4 Adenylating 3' Ends\*



## **RNAseq: Library Preparation**

#### Figure 5 Ligating Adapters





## **RNAseq: Library Preparation**

Figure 6 Enriching DNA Fragments\*






Figure 7 LS Final Library





### RNA Library QC: Electrophoresis and Electropherogram

### **Electrophoresis**




### Electropherogram







## RNAseq: Library QC

Libraries are evaluated using a bioanalyzer or a tape station to create an electropherogram to assess quality







When assessing library quality look for the following:


- NO adapter dimers! Why?
- Library size is consistent with the number of desired sequencing cycles
  - ➤ If you're sequencing at PE 100, what is a good library size? (hint: library size = insert length + adapter length)


Assuming adapters are ~65bp each (130bp): ~300bp that would give an insert length of ~170, allowing a ~30bp overlap between R1 and R2


What is the size range of library 1? Library 2? ~175-250bp; ~150-275bp, Modified from UC Davis Bioinformatics RNAseq Training

## RNAseq: Library QC

- Assess library quality (bioanalyzer, TapeStation)
  - Adapter dimers
  - > Fragment size
- Determine library quantity (Qubit, qPCR)







% AD	% PF	% AD Reads		
Control	69.54	0.24		
10%	10.87	84.25		
5%	21.39	60.44		
1%	51.88	6.46		



### **Cluster Formation**






Figure 4. Bridge PCR - a PCR method used to amplify samples for sequencing.

## Illumina's High Throughput Sequencing by Synthesis



- For paired-end (PE) sequencing, after read
   1 is sequenced, forward strand reagents
   are washed
- The index read(s) are sequenced next
- Sequencing of the reverse read (read 2) is initiated after the index read(s)

Figure 9 Read 2 Sequencing



# **Sequencing Parameters**

	Read length	Sequencing depth	Paired-end (PE) or Single- end (SE)
DNAseq	<ul> <li>Longer reads enable greater confidence in taxonomic classifications and functional annotations</li> <li>"Assembly" is often performed with short reads to facilitate this</li> <li>GL standard is 2x250bp</li> </ul>	<ul> <li>Greater depth increases the likelihood of sequencing low-abundance organisms (if metagenomics) and detecting things like single-nucleotide variants and genetic rearrangements with greater confidence</li> <li>GL standards:         <ul> <li>Single organism or tissue:</li> <li>Re-sequencing (reference available): 10X minimum</li> <li>De novo sequencing (no reference available): 50X minimum</li> <li>Metagenomics (mixed community): 10M per sample, minimum</li> </ul> </li> </ul>	PE is generally preferred
RNAseq	<ul> <li>Longer reads increase gene ID confidence</li> <li>GL standard is 2x150bp for bulk RNAseq</li> </ul>	<ul> <li>Greater depth increases the likelihood of sequencing low-abundant transcripts, detecting novel transcripts, and quantifying isoforms</li> <li>Greater depth is necessary for ribo-depleted samples (vs. poly-A enriched samples) – for RNAseq</li> <li>GL RNAseq standard for mammalian samples prepared with ripo-depletion is ~60M reads/sample and ~40M reads/sample for samples prepared with polyA-selection.</li> <li>More replicates is usually preferred over greater depth</li> </ul>	➤ PE is preferred

## Illumina's High Throughput Sequencing by Synthesis

https://www.youtube.com/watch?v=fCd6B5HRaZ8



